

# bLINK


## Integration Instructions for Edge Cloud



*Edge Cloud's integration allows bLINK to sync with your schedule so bLINK can send appointment reminders.*

## Setup

To enable the bLINK integration, open the Editors from the Tools tab of the Home ribbon bar. Click System Integrations under the Integrations heading.



Under the bLINK by Solutions by Design heading, click the **Enable** button.

## Contact

**bLINK**  
Solutions by Design  
451 Clovis Avenue, Suite 130  
Clovis, CA 93612  
[www.solutionsbydesign.com/blink-pro.php](http://www.solutionsbydesign.com/blink-pro.php)  
800.888.4084

**Edge Cloud**  
Ortho2  
1107 Buckeye Avenue  
Ames, IA 50010  
[www.ortho2.com](http://www.ortho2.com)  
800.346.4504 or 515.233.1026